

Recombinant Human EIF1AX protein ab101058

1 图像

| | |
|--|---|
| 描述 | |
| 产品名称 | 重组人EIF1AX蛋白 |
| 纯度 | > 90 % SDS-PAGE. ab101058 was purified by using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA. |
| 表达系统 | Escherichia coli |
| Accession | <u>P47813</u> |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Human |
| 序列 | MGSSHHHHHH SSGLVPRGSH MPKNKGKGGK NRRRGKNENE SEKRELVFKE DGQEYAQVIK MLGNGRLEAM CFDGVKRLCH IRGKLKKVW INTSDIILVG LRDYQDNKAD VILKYNADDA RSLKAYGELP EHAKINETDT FGPDDDEIQ FDDIGDDDED IDDI |
| 预测分子量 | 19 kDa including tags |
| 氨基酸 | 1 to 144 |
| 标签 | His tag N-Terminus |
| 技术指标 | |
| Our Abpromise guarantee covers the use of ab101058 in the following tested applications. | |
| The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user. | |
| 应用 | Mass Spectrometry SDS-PAGE |
| 质谱法 | MALDI-TOF |
| 形式 | Liquid |
| 制备和贮存 | |

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.077% DTT, 0.316% Tris HCl, 40% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

常规信息

相关性

This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA_i and eIF3) to the 5' end of capped RNA.

图片



15% SDS-PAGE showing ab101058 (3 µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors