

Recombinant Human Elafin/ESI protein ab153768

描述

产品名称	重组人Elafin/ESI蛋白
纯度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by reducing SDS-PAGE. ab153768 is supplied as a 0.2 µm filtered solution.
内毒素水平	< 1.000 Eu/µg
表达系统	HEK 293 cells
Accession	<u>P19957</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	AVTGVPVKGQDTVKGRVPFNGQDPVKGQVSVKGQDKVKAQEP VKGPVSTK PGSCPIILIRCAMLNPPNRCLKDTDCPGIKKCEGSCGMACF VPQ
预测分子量	10 kDa
氨基酸	23 to 117
标签	His tag C-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab153768** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	Previously labelled as Elafin/Skalp.

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50
--------	---

Constituents: 0.32% Tris HCl, 0.88% Sodium chloride

常规信息

功能	Neutrophil and pancreatic elastase-specific inhibitor of skin. It may prevent elastase-mediated tissue proteolysis.
序列相似性	Contains 1 WAP domain.
结构域	Consists of two domains: the transglutaminase substrate domain (cementin moiety) and the elastase inhibitor domain. The transglutaminase substrate domain serves as an anchor to localize elafin covalently to specific sites on extracellular matrix proteins.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors