

Recombinant Human Eotaxin 2 protein ab9834

描述

产品名称	重组人Eotaxin 2蛋白
生物活性	Determined by its ability to chemoattract using human peripheral blood eosinophils using a concentration of 50.0-100.0 ng/ml.
纯度	> 98 % SDS-PAGE. Sterile filtered Greater than 98% pure by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).
内毒素水平	< 0.100 Eu/μg
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	VVIPSPCCMF FVSKRIPENR VVSYQLSSRS TCLKGGVIFT TTKGQQFCGD PKQEWVQRYM KNLDAKQKKA SPRARAVA
预测分子量	13 kDa

技术指标

Our **Abpromise guarantee** covers the use of **ab9834** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

形式	Lyophilized
补充说明	The biological activity of this product is determined by its ability to chemoattract human peripheral blood eosinophils using a concentration range of 50.0-100.0 ng/ml.

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
复溶	For lot specific reconstitution information please contact our Scientific Support Team.

常规信息

功能	Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.
组织特异性	Activated monocytes and activated T lymphocytes.
序列相似性	Belongs to the intercrine beta (chemokine CC) family.
翻译后修饰	N-glycosylated.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors