


Recombinant Human FABP4 protein ab63283

2 图像


描述	
产品名称	重组人FABP4蛋白
纯度	> 90 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in cell culture.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	Yes
性质	Recombinant
种属	Human
序列	M-CDAFVGTWK LVSSNFDDY MKEVGVGFAT RKVAGMAKPN MIISVNGDVI TIKSESTFKN TEISFILGQE FDEVTADDRK VKSTITLDGG VLVHVQKWDG KSTTIKREKRE DDKLVVECVM KGV TSTRVYE RA
技术指标	
Our Abpromise guarantee covers the use of ab63283 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot ELISA
形式	Lyophilized
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituent: PBS
复溶	Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.
常规信息	

功能	Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.
序列相似性	Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.
结构域	Forms a beta-barrel structure that accommodates hydrophobic ligands in its interior.
细胞定位	Cytoplasm. Nucleus. Depending on the nature of the ligand, a conformation change exposes a nuclear localization motif and the protein is transported into the nucleus. Subject to constitutive nuclear export.

图片



LANE 1. M.W. marker - 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane.



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors