

Recombinant Human Fc epsilon RI/FCER1A protein ab114334

1 图像

描述	
产品名称	重组人Fc epsilon RI/FCER1A蛋白
表达系统	Wheat germ
Accession	<b>P12319</b>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAPAMESPTLLCVALLFFAPDGVLAVPQKPKVSLNPPWNRIF KGENVTLT CNGNNFFEVSSTKWFFHNGSLSEETNSSLNIVNAKFEDSGEYK CQHQQVNE SEPVYLEVFSDWLLLQASAEVVMEGQPLFLRCHGWRNWDVYK VIYYKDGE ALKYWYENHNISITNATVEDSGTYCYCTGKVVQLDYESEPLNI TVIKAPRE KYWLQFFIPLLVVILFAVDTGLFISTQQQVTFLLKIKRTRKG FRLLNPHP KPNPKNN
预测分子量	54 kDa including tags
氨基酸	1 to 257
标签	GST tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab114334</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	ELISA SDS-PAGE Western blot
形式	Liquid
补充说明	This product was previously labelled as Fc epsilon RI.

## 制备和贮存

### 稳定性和存储

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

## 常规信息

### 功能

Binds to the Fc region of immunoglobulins epsilon. High affinity receptor. Responsible for initiating the allergic response. Binding of allergen to receptor-bound IgE leads to cell activation and the release of mediators (such as histamine) responsible for the manifestations of allergy. The same receptor also induces the secretion of important lymphokines.


### 序列相似性

Contains 2 Ig-like (immunoglobulin-like) domains.

### 细胞定位

Cell membrane.

## 图片



ab114334 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors