

Recombinant human FGF18 protein ab50129

描述	
产品名称	重组人FGF18蛋白
生物活性	Biological Activity: Determined by dose-dependent ability to reduce tetrazolium salt, WST-8, by dehydrogenase activities of BaF3 cells expressing FGF receptors using Cell Counting Kit-8 (CCK-8).
纯度	> 95 % SDS-PAGE. Greater than 95% by SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MAEENVDFRI HVENQTRARD DVSRKQLRLY QLYSRTSGKH IQVLGRRISA RGEDGDKYAAQ LLVETDTFGS QVRIKGKETE FYLCMNRKGK LVGKPDGTSK ECVFIEKVLE NNYTALMSAK YSGWYVGFTK KGRPRKGPKT RENQQDVHFM KRYPKGQPEL QKPFKYTTVT KRSRRIRPTH PA
技术指标	
Our Abpromise guarantee covers the use of ab50129 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Functional Studies
形式	Lyophilized
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituents: 0.121% Tris, 0.116% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.

复溶	Centrifuge the vial prior to opening. Reconstitute in 5mM Tris, pH 8.0 to a concentration of 0.1-1.0 mg/ml. Do not vortex. This solution can be stored at 2-8°C; for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20°C; to -80°C;.
----	---

常规信息

功能	Plays an important role in the regulation of cell proliferation, cell differentiation and cell migration. Required for normal ossification and bone development. Stimulates hepatic and intestinal proliferation.
序列相似性	Belongs to the heparin-binding growth factors family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors