

Recombinant human FKBP25 protein ab93683

1 图像

描述		
产品名称	重组人FKBP25蛋白	
生物活性	Biological activity: Specific activity is > 900nmol/min/mg, and is defined as the amount of enzyme that cleaves 1nmol of suc-AAPF-pNA per minute at 37C in Tris-HCl pH 8.0 using chymotrypsin	
纯度	> 90 % SDS-PAGE. ab93683 is purified using conventional chromatography techniques.	
表达系统	Escherichia coli	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MAAAVPQRAWTV EQLRSEQLPKKDIIKFLQEHGSDSFLAEHK LLGNIKNV AKTANKDHLVTAYNHLFETKRFGTESISKVSEQVKNVKLNE DKPKETKS EETLDEGPPKYTKSVLKKGDKTNFPKKGDVVHCWYTGT LQDG TVFDTNIQ TSAKKKKNAKPSFKVGVGKVIRGWDEALLTMSKGEKARLEIE PEWAYGKK GQPDAKIPPNAKLTFEVELVDID	

技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab93683 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Functional Studies SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

pH: 8.00


Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

This product is an active protein and may elicit a biological response in vivo, handle with caution.

常规信息

功能	FK506- and rapamycin-binding proteins (FKBPs) constitute a family of receptors for the two immunosuppressants which inhibit T-cell proliferation by arresting two distinct cytoplasmic signal transmission pathways. PPlases accelerate the folding of proteins.
序列相似性	Belongs to the FKBP-type PPlase family. Contains 1 PPlase FKBP-type domain.
细胞定位	Nucleus.

图片



15% SDS-PAGE showing ab93683 at approximately 25 kDa (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors