

Recombinant Human FLR protein ab79183

1 图像

描述	
产品名称	重组人FLR蛋白
纯度	> 95 % SDS-PAGE. ab79183 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAVKKIAIFG ATGQTGLTTL AQAQAGYEV TVLVRDSSRL PSEGPRPAHV VVGDLQAAD VDKTVAGQDA VIVLLGTRND LSPTTVMSEG ARNIVAAMKA HGVDKVVACT SAFLLWDPTK VPPRLQAVTD DHIRMHKVLRL ESGLKYVAVM PPHIGDQPLT GAYTVTLQDGR GPSRVISKHD LGHFMLRCLT TDEYDGHSTY PSHQYQ


技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab79183</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as BLVRB

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  pH: 8.50 Constituents: 0.0154% DTT, 0.242% Tris, 10% Glycerol (glycerin, glycerine)

## 常规信息

功能	Broad specificity oxidoreductase that catalyzes the NADPH-dependent reduction of a variety of flavins, such as riboflavin, FAD or FMN, biliverdins, methemoglobin and PQQ (pyrroloquinoline quinone). Contributes to heme catabolism and metabolizes linear tetrapyrroles. Can also reduce the complexed Fe(3+) iron to Fe(2+) in the presence of FMN and NADPH. In the liver, converts biliverdin to bilirubin.
组织特异性	Predominantly expressed in liver and erythrocytes. At lower levels in heart, lung, adrenal gland and cerebrum.
细胞定位	Cytoplasm.

## 图片



15% SDS-PAGE showing ab79183 at approximately 22kDa (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors