

## Recombinant Human GBP1 protein ab114960

1 图像

描述

<b>产品名称</b>	重组人GBP1蛋白
<b>表达系统</b>	Wheat germ
<b>Accession</b>	<b><u>P32455</u></b>
<b>蛋白长度</b>	Full length protein
<b>无动物成分</b>	No
<b>性质</b>	Recombinant
<b>种属</b>	Human
<b>序列</b>	<pre> MASEIHMTGPMCLIENTNGRLMANPEALKILSAITQPMVVVA IVGLYRTG KSYLMNKLAGKKKGFSLGSTVQSHTKGIWMWCVPHPKKPGHI LVLLDTEG LGDVEKGDNDQNSWIFALAVLLSSTFVYNSIGTINQQAMDQL YYVTELTH RIRSKSSPDENENEVEDSADFVSFFPDFVWTLRDFSLDLEAD GQPLTPDE YLTYSCLKKKGTSQKDETFNLPRLCIRKFFPKKKCFVDRPV HRRKLAQL EKLQDEELDPEFVQQVADFCSYIFSNSKTKTSLGGIQVNGPR LESLVLTY VNAISSGDLPCMENAVLALAQIENSAAVQKAI AHYEQQMGQK VQLPTESL QELLDLHRDSEREAIEVFIRSSFKDVDHLFQKELAAQLEKKR DDFCKQNQ EASSDRCSGLLQVIFSPLEEEVKAGIYSKPGGYRLFVQKLQD LKKKYEE PRKGIQAEEILQTYLKSKE SMTDAILQTDQTLTEKEKEIEVE RVKAES AQ ASAKMLQEMQRKNEQMMEQKERSYQEHLKQLTEKMENDRVQL LKEQERTL ALKLQEQEQLLKEGFQKESRIMKNEIQDLQTKMRRRRACTIS </pre>
<b>预测分子量</b>	91 kDa including tags
<b>氨基酸</b>	1 to 592

## 这个指标

---

Our **Abpromise guarantee** covers the use of **ab114960** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	SDS-PAGE
	Western blot
形式	Liquid

## 制备和贮存

---

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

## 常规信息


---

功能	Binds GTP, GDP and GMP.
序列相似性	Belongs to the GBP family.
细胞定位	Cell membrane.

---

## 图片

---



12.5% SDS-PAGE showing ab114960 at approximately 91.23kDa.  
Stained with Coomassie Blue.

SDS-PAGE - Recombinant Human GBP1 protein  
(ab114960)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors