

Product datasheet

Recombinant Human GCIP interacting protein p29 ab113599

1 图像

描述

产品名称	重组人GCIP interacting蛋白p29
纯度	> 90 % SDS-PAGE. ab113599 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.
表达系统	Escherichia coli
Accession	<u>O95926</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMGSDYEKVKLLEISAEDAER WERKKKRKNP DLGFSDYAAAQLRQYHRLTKQIKPDMETYERLREKHGEEFFP TSNSLLHG THVPSTEEIDRMVIDLEKQIEKRDKYSRRRPYNDDADIDYIN ERNAKFNK KAERFYGKYTAEIKQNLERGTA
预测分子量	21 kDa including tags
氨基酸	94 to 243
标签	His tag N-Terminus
描述	重组人GCIP interacting蛋白p29

技术指标

Our **Abpromise guarantee** covers the use of **ab113599** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid


稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
 pH: 8.00
 Constituents: 0.08% DTT, 0.32% Tris HCl, 0.06% EDTA, 50% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

常规信息

功能 May be involved in pre-mRNA splicing.
组织特异性 Abundantly expressed in the heart, skeletal muscle and kidney. Expressed at lower level other tissues.
序列相似性 Belongs to the SYF2 family.
细胞定位 Nucleus.

图片



15% SDS-PAGE showing ab113599 at approximately 20.5 kDa (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors