

Recombinant Human GNG11 protein (His tag) ab177606

1 图像

描述		
产品名称	重组人GNG11蛋白(His tag)	
纯度	> 90 % SDS-PAGE. ab177606 was purified using conventional chromatography techniques.	
表达系统	Escherichia coli	
Accession	<u>NP_004117</u>	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MGSSHHHHHHSSGLVPRGSHMGSM PALHIEDLPEKEKLKMEV EQLRKEVK LQRQQVSKCSEEEKNYIEERSGEDPLVKGIPEDKNPFKEKGS C	
预测分子量	11 kDa including tags	
氨基酸	1 to 70	
标签	His tag N-Terminus	
额外的序列信息	(NP_004117)	

技术指标		
Our Abpromise guarantee covers the use of ab177606 in the following tested applications.		
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.		
应用	Mass Spectrometry SDS-PAGE	
质谱法	MALDI-TOF	
形式	Liquid	

制备和贮存		
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -	

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 10% Glycerol (glycerin, glycerine), 0.02% DTT

常规信息

功能	Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction.
组织特异性	Abundantly expressed in all tissues tested except for brain.
序列相似性	Belongs to the G protein gamma family.
细胞定位	Cell membrane.

图片



15% SDS-PAGE analysis of ab177606 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors