

Recombinant Human GSTA1 protein ab95381

2 图像

描述	
产品名称	重组人GSTA1蛋白
纯度	> 90 % SDS-PAGE. purified by using conventional chromatography techniques
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAEKPKLHYF NARGRMESTR WLLAAAGVEF EEKFIKSAED LDKLRNDGYL MFQQVPMVEI DGMKLVQTRA ILNYIASKYN LYGKDIKERA LIDMYIEGIA DLGEMILLLP VCPPEEKDAK LALIKEKIKN RYFPAFEKVL KSHGQDYLVG NKLSRADIHL VELLYYVEEL DSSLISSFPL LKALKTRISN LPTVKKFLQP GSPRKPPMDE KSLEEARKIF RF

技术指标	
Our Abpromise guarantee covers the use of ab95381 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Western blot
形式	Liquid
补充说明	Previously labelled as Glutathione S Transferase alpha 1. Previously labelled as Glutathione S Transferase alpha 1.

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

常规信息


功能	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
组织特异性	Liver.
序列相似性	Belongs to the GST superfamily. Alpha family. Contains 1 GST C-terminal domain. Contains 1 GST N-terminal domain.
结构域	The C-terminal domain may form a component of the hydrophobic substrate-binding site, but in contrast appears not to be directly involved in GSH binding and is not absolutely essential for catalytic activity.
细胞定位	Cytoplasm.

图片



15% SDS-PAGE image (3ug) generated using ab95381

SDS-PAGE - Recombinant Human GSTA1 protein
(ab95381)



Western blot - Recombinant Human GSTA1 protein
(ab95381)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors