



Product datasheet

Recombinant Human HIGD1A protein ab125988

描述	
产品名称	重组人HIGD1A蛋白
纯度	> 95 % SDS-PAGE. Purified via His tag
表达系统	Escherichia coli
Accession	<u>Q9Y241</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MSTDTGVSLPSYEEDQGSKLIRKAKEAPFVPVGIAGFAAIVA YGLYKLKS RGNTKMSIHLIHMRVAAQGFVVGAMTVGMGYSMYREFWAKPK P
氨基酸	1 to 93

技术指标	
Our Abpromise guarantee covers the use of ab125988 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Lyophilized

制备和贮存	
稳定性和存储	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
复溶	Reconstitute with water to desired concentration.

常规信息	
序列相似性	Contains 1 HIG1 domain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors