

Recombinant Human HLF protein ab113145

1 图像

描述	
产品名称	重组人HLF蛋白
纯度	> 85 % SDS-PAGE. > 85% by SDS - PAGE. ab113145 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	<u>Q16534</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMEKMSRPLPLNPTFIPPPYGVL RSLLLENPL KLPLHHEDAFSKDKDKEK KLDDESNSPTVPQSAFLGPTLWDKTLPYDG DTFQLEYMDLEEFLSENGIPPSPSQHDHSPHPPGLQPASSAA PSV MDL SSRASAPLHPGIPSPNCMQSPIRPGQLLPANRNTPSPIDPDT IQVPVGYE PDPADLALSSIPGQEMFDPRKR KFSEEELKPQPMIKKARKVFIPDDLK DDKYWARRRKNMAAKRSRDARRLKENQIAIRASFLEKENSA LRQEVAD LRKELGKCKNILAKYEARHGPL
预测分子量	35 kDa including tags
氨基酸	1 to 295
标签	His tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab113145</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid

## 制备和贮存

### 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.002% PMSF, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

## 常规信息

### 组织特异性

Highly expressed in liver; lower levels in lung and kidney.

### 疾病相关

Note=A chromosomal aberration involving HLF is a cause of pre-B-cell acute lymphoblastic leukemia (B-ALL). Translocation t(17;19)(q22;p13.3) with TCF3.

### 序列相似性


Belongs to the bZIP family. PAR subfamily.

Contains 1 bZIP domain.

### 细胞定位

Nucleus.

## 图片



15% SDS-PAGE analysis of ab113145 (3µg)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors