

Recombinant Human Hsp40 protein ab51305

1 References 1 图像

描述	
产品名称	重组人Hsp40蛋白
纯度	> 95 % SDS-PAGE. HSP40 (amino acids 1-376) was overexpressed in <i>E. coli</i> and purified to apparent homogeneity by using conventional column chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAKQDYIEIL GVSKTAEHE IRKAYKRLAM KYHPDRNQGD KEAEAKFKEI KEAYEVLTD QKRAAYDQYG HAAFEQGGMG GGGFGGGADF SDIFGDFVGD IFGGGRGRQR AARGADLRN MELTLEEAVR GVTKEIRIPT LEECDVCHGS GAKPGTQPQT CPTCHGSGQV QMRQGFFAVQ QTCPHCQGRG TLIKDPCNKC HGHGRVERSK TLSVKIPAGV DTGDRIRLAG EGEAGEHGAP AGDLYVQVQV KQHPIFEREG NNLYCEVPIN FAMAALGGEI EVPTLDGRVK LKVPGETQTG KLFRMRGKGV KSVRGAQGD LLCRVVVETP VGLNERQKQL LQELQESFGG PTGEHNSPRS KSFFDGVKKF FDDLTR
氨基酸	1 to 376
技术指标	
Our Abpromise guarantee covers the use of ab51305 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid
制备和贮存	

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

pH: 8.8

Constituents: 0.077% DTT, 0.395% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

功能

Interacts with HSP70 and can stimulate its ATPase activity. Stimulates the association between HSC70 and HIP.

序列相似性

Contains 1 J domain.

细胞定位

Cytoplasm. Nucleus. Nucleus > nucleolus. Translocates rapidly from the cytoplasm to the nucleus, and especially to the nucleoli, upon heat shock.

图片



ab51305 in 13% SDS-PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors