



Product datasheet

Recombinant Human IL-17B protein ab50161

|  |  |
|--|--|
| 描述   |  |
| 产品名称   | 重组人IL-17B蛋白  |
| 纯度   | > 98 % SDS-PAGE.<br>Greater than 98% by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).  |
| 表达系统   | Escherichia coli   |
| 蛋白长度   | Full length protein  |
| 无动物成分  | No   |
| 性质   | Recombinant  |
| 种属   | Human  |
| 序列   | MQPRSPKSKR KGQGRPGPLA PGPHQVPLDL<br>VSRMKPYARM EEYERNIEEM VAQLRNSSEL<br>AQRKCEVNLQ LWMSNKRSLs PWGYSINHDP<br>SRIPVDLPEA RCLCLGCVNP FTMQEDRSMV<br>SVPVFSQVPV RRRLCPPPPR TGPCRQRAVM<br>ETIAVGCTCI F |
| 技术指标   |  |
| Our <b>Abpromise guarantee</b> covers the use of <b>ab50161</b> in the following tested applications.                                |  |
| The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user. |  |
| 应用   | SDS-PAGE<br>Functional Studies   |
| 形式   | Lyophilized  |
| 制备和贮存  |  |
| 稳定性和存储   | Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.   |
| 复溶   | For lot specific reconstitution information please contact our Scientific Support Team.  |
| 常规信息   |  |

|       |   |
|-------|---|
| 功能    | Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.   |
| 组织特异性 | Expressed in adult pancreas, small intestine, stomach, spinal cord and testis. Less pronounced expression in prostate, colon mucosal lining, and ovary. |
| 序列相似性 | Belongs to the IL-17 family.  |
| 细胞定位  | Secreted.   |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors