

Recombinant Human Kallikrein 13 protein (denatured) ab202250

1 图像

描述	
产品名称	重组人Kallikrein 13蛋白(denatured)
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>Q9UKR3</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MGSGGVSQES SKVLNTNGTS GFLPGGYTCF PHSQPWQAAL LVQGRLLCGG VLVHPKWVLT AAHCLKEGLK VYLGKHALGR VEAGEQVREV VHSIPHPEYR RSPTHLNHDH DIMLLELQSP VQLTGYIQL PLSHNNRLTP GTTCRVSGWG TTTSPQVNYP KTLQCANIQL RSDEECRQVY PGKITDNMLC AGTKEGGKDS CEGDSGGPLV CNRTLYGIVS WGDFPCGQPD RPGVYTRVSR YVLWIRETIR KYETQQQKWL KGPQ
预测分子量	31 kDa including tags
氨基酸	17 to 277
标签	His tag N-Terminus
额外的序列信息	NP_056411. This product is for the mature full length protein from aa 17 to 277. The signal peptide is not included.

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab202250</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid

## 制备和贮存

### 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 10% Glycerol (glycerin, glycerine), 0.32% Tris HCl

## 常规信息

### 组织特异性

Expressed in prostate, breast, testis and salivary gland.

### 序列相似性


Belongs to the peptidase S1 family. Kallikrein subfamily.

Contains 1 peptidase S1 domain.

### 细胞定位

Secreted.

## 图片



15% SDS-PAGE analysis of ab202250 (3 µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors