

Recombinant Human KDM4A / JHDM3A / JMJD2A protein ab125541

1 图像

描述

产品名称	重组人KDM4A / JHDM3A / JMJD2A蛋白
纯度	> 90 % Densitometry. Purity was determined to be >90% by densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<u>O75164</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	150 kDa including tags
氨基酸	1 to 886
标签	GST tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab125541** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
--------	---

常规信息

功能	Histone demethylase that specifically demethylates 'Lys-9' and 'Lys-36' residues of histone H3, thereby playing a central role in histone code. Does not demethylate histone H3 'Lys-4', H3 'Lys-27' nor H4 'Lys-20'. Demethylates trimethylated H3 'Lys-9' and H3 'Lys-36' residue, while it has no activity on mono- and dimethylated residues. Demethylation of Lys residue generates formaldehyde and succinate. Participates in transcriptional repression of ASCL2 and E2F-responsive promoters via the recruitment of histone deacetylases and NCOR1, respectively. Isoform 2: Crucial for muscle differentiation, promotes transcriptional activation of the Myog gene by directing the removal of repressive chromatin marks at its promoter. Lacks the N-terminal demethylase domain.
组织特异性	Ubiquitous.
序列相似性	Belongs to the JHDM3 histone demethylase family. Contains 1 C2HC pre-PHD-type zinc finger. Contains 1 JmjC domain. Contains 1 JmjN domain. Contains 2 PHD-type zinc fingers. Contains 2 Tudor domains.
结构域	The 2 Tudor domains recognize and bind methylated histone H3 'Lys-4' residue (H3K4me). Double Tudor domain has an interdigitated structure and the unusual fold is required for its ability to bind methylated histone tails. Trimethylated H3 'Lys-4' (H3K4me3) is bound in a cage of 3 aromatic residues, 2 of which are from the Tudor domain 2, while the binding specificity is determined by side-chain interactions involving residues from the Tudor domain 1. The Tudor domains are also able to bind trimethylated histone H3 'Lys-9' (H3K9me3), di- and trimethylated H4 'Lys-20' (H4K20me2 and H4K20me3). Has high affinity for H4K20me2, blocking recruitment of proteins such as TP53BP1.
翻译后修饰	Ubiquitinated by RNF8 and RNF168 following DNA damage, leading to its degradation. Degradation promotes accessibility of H4K20me2 mark for DNA repair protein TP53BP1, which is then recruited.
细胞定位	Nucleus.

图片



SDS-PAGE analysis of ab125541.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors