

# Recombinant Human LCN1 protein ab116833

[1 References](#) [1 图像](#)

### 描述

产品名称	重组人LCN1蛋白
表达系统	Wheat germ
Accession	<b><u>P31025</u></b>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MKPLLLAVSLGLIAALQAHLLASDEEIQDVSGTWYLKAMTV DREFPEMN LESVTPMTLTTLEGGNLEAKVTMLISGRCQEVKAVLEKTDEP GKYTADGG KHVAYIIRSHVKDHYIFYCEGELHGKPVRGVKLVGRDPKNNL EALEDFEK AAGARGLSTESILIPRQSETCSPGSD
预测分子量	45 kDa including tags
氨基酸	1 to 176

### 技术指标

Our **Abpromise guarantee** covers the use of **ab116833** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid

### 制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

## 常规信息


---

<b>功能</b>	Could play a role in taste reception. Could be necessary for the concentration and delivery of sapid molecules in the gustatory system. Can bind various ligands, with chemical structures ranging from lipids and retinoids to the macrocyclic antibiotic rifampicin and even to microbial siderophores. Exhibits an extremely wide ligand pocket.
<b>组织特异性</b>	Mainly expressed in lachrymal and salivary glands. Also expressed in the prostate.
<b>序列相似性</b>	Belongs to the calycin superfamily. Lipocalin family.
<b>细胞定位</b>	Secreted.

---

## 图片

---



12.5% SDS-PAGE stained with Coomassie Blue showing ab116833 at approximately 45.43 kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors