

Recombinant Human LSM1 protein ab101157

1 图像

描述	
产品名称	重组人LSM1蛋白
纯度	> 95 % SDS-PAGE. ab101157 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>O15116</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMNYMPGTASLIEDIDKKHLVLL RDGRTLIG FLRSIDQFANLVLHQTVRIHVGKKYGDIPRGIFVVRGENVV LLGEIDLE KESDTPLQQVSIEEILEEQRVEQQTKEAEKLVQALKDRGL SIPRADTL DEY
预测分子量	17 kDa including tags
氨基酸	1 to 133
标签	His tag N-Terminus
技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab101157 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.


pH: 7.50

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

功能	Plays a role in replication-dependent histone mRNA degradation. Binds specifically to the 3'-terminal U-tract of U6 snRNA.
组织特异性	Has elevated expression in pancreatic cancer and in several cancer-derived cell lines.
序列相似性	Belongs to the snRNP Sm proteins family.
细胞定位	Nucleus.

图片



15% SDS-PAGE showing ab101157 at approximately 17.3kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors