

Recombinant Human LSM4 protein ab101871

1 图像

描述

产品名称	重组人LSM4蛋白
纯度	> 90 % SDS-PAGE. ab101871 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	Q9Y4Z0
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMLPLSLLKTAQNHPMLVELK NGETYNGHLV SCDNWMNINLREVICTSRDGDKFWRMPECYIRGSTIKYLRI DEIIDMVK EEVVAKGRGRGGLQQKQQKGRGMGGAGRGVFGGRGRGGIPG TGRGQPEK KPGRQAGKQ
预测分子量	18 kDa including tags
氨基酸	1 to 139
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab101871** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

功能	Binds specifically to the 3'-terminal U-tract of U6 snRNA.
序列相似性	Belongs to the snRNP Sm proteins family.
细胞定位	Nucleus.

图片



15% SDS PAGE analysis of 3µg ab101871.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors