

Recombinant Human MED4 protein ab95264

1 图像

描述	
产品名称	重组人MED4蛋白
纯度	> 80 % SDS-PAGE. Purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAASSSGEKE KERLGGLGV AGGNSTRERL LSALEDLEVL SRELIEMLAISRNRKLLQAG EENQVLELLIHRDGEFQELMKLALNQGKIH HEMQVLEKEVEKRDGDIQQLQKQLKEAEQI LATAVYQAKEKLKSIEKARKGAISSEEEIK YAHRIASNAVCAPLTWVPGDPRRPYPTDL EMRSGLLGQMNNPSTNGVNGHLPGDALAAG RLPDVLAPQYPWQSNDSMNMLPPNHSSDF LLEPPGHNKEDEDDVEIMSTDSSSSSESDELEHHHHHH

技术指标	
Our Abpromise guarantee covers the use of ab95264 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
形式	Liquid


制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium

chloride

常规信息

功能	Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors.
序列相似性	Belongs to the Mediator complex subunit 4 family.
细胞定位	Nucleus.

图片



15% SDS-PAGE analysis of 3µg ab95264.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors