

Recombinant Human MIP-3 beta/CCL19 protein ab97899

1 图像

描述	
产品名称	重组人MIP-3 beta/CCL19蛋白
纯度	> 95 % SDS-PAGE. ab97899 is purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.
内毒素水平	< 1.000 Eu/μg
表达系统	Escherichia coli
Accession	<u>Q99731</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<i>MASMTGGQQM GRGSHMG</i> TND AEDCCLSVTQ KPIPGYIVRN FHYLLIKDGC RVPVVFTTL RGRQLCAPPD QPWVERIIQR LQRTSAKMKR RSS
预测分子量	10 kDa including tags
氨基酸	22 to 98
标签	T7 tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab97899 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF-TOF
形式	Liquid
补充说明	This product was previously labelled as Macrophage Inflammatory Protein 3 beta

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 7.40

Constituents: PBS, 10% Glycerol (glycerin, glycerine)

常规信息

功能

May play a role not only in inflammatory and immunological responses but also in normal lymphocyte recirculation and homing. May play an important role in trafficking of T-cells in thymus, and T-cell and B-cell migration to secondary lymphoid organs. Specifically binds to chemokine receptor CCR7. Recombinant CCL19 shows potent chemotactic activity for T-cells and B-cells but not for granulocytes and monocytes.

组织特异性

Expressed at high levels in the lymph nodes, thymus and appendix. Intermediate levels seen in colon and trachea, while low levels found in spleen, small intestine, lung, kidney and stomach.


序列相似性

Belongs to the intercrine beta (chemokine CC) family.

细胞定位

Secreted.

图片



15% SDS-PAGE showing ab97899 at approximately 10.4kDa (3μg).

SDS-PAGE - Recombinant Human MIP-3
beta/CCL19 protein (ab97899)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors