

Recombinant human MST4 protein ab60341

4 图像

描述

产品名称	重组人MST4蛋白
纯度	> 90 % SDS-PAGE. Purity: >90% as determined by densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human

技术指标

Our **Abpromise guarantee** covers the use of **ab60341** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE Western blot
形式	Liquid
补充说明	ab43614 (Human Myelin Basic Protein full length protein) can be utilized as a substrate for assessing Kinase activity


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.0038% EGTA, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

常规信息


功能	Mediator of cell growth. Modulates apoptosis.
序列相似性	Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. STE20 subfamily. Contains 1 protein kinase domain.
细胞定位	Cytoplasm. Golgi apparatus membrane. Localizes to the Golgi apparatus.

图片




The specific activity of MST4 (ab60341) was determined to be 23 nmol/min/mg as per activity assay protocol

Functional Studies - Recombinant human MST4 protein (ab60341)



SDS PAGE analysis of ab60341

SDS-PAGE - Recombinant human MST4 protein (ab60341)



Sample Kinase Activity Plot.

Functional Studies - Recombinant human MST4 protein (ab60341)



Western blot - Recombinant human MST4 protein (ab60341)

Anti-MST4 antibody [EP1864Y] (**ab52491**) at 1/20000 dilution + Recombinant human MST4 protein (ab60341) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (**ab97080**) at 1/5000 dilution

Exposure time: 20 seconds

Ab52491 recognizes the full length tagged recombinant MST4 protein (ab60341) which has an expected molecular weight of 72 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors