

Recombinant Human Myostatin Propeptide protein ab63300

1 图像

描述	
产品名称	重组人Myostatin Propeptide蛋白
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>O14793</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MGNENSEQKE NVEKEGLCNA CTWRQNTKSS RIEAIKIQIL SKLRLETAPN ISKDVI RQLL PKAPPLRELI DQYDVQRDDS SDGSLEDDDY HATTETIITM PTESDFLMQV DGKPKCCFFK FSSKIQYNKV VKAQLWIYLR PVETPTTVFV QILRLIKPMK DGTRYTGIRS LKLDMNPGTG IWQSIDVKTV LQNLWKQPES NLGIEIKALD ENGHDLAVTF PGPGEDGLNP FLEVKVTDTP KRSRRKLN
氨基酸	24 to 266

技术指标	
Our Abpromise guarantee covers the use of ab63300 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot ELISA
形式	Lyophilized
补充说明	Product is NOT sterile! Please filter the product by an appropriate sterile filter before using it in cell culture.

制备和贮存	
-------	--

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
复溶	Add deionized water to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely.
常规信息	
功能	Acts specifically as a negative regulator of skeletal muscle growth.
疾病相关	Defects in MSTN are the cause of muscle hypertrophy (MSLHP) [MIM:614160]. MSLHP is a condition characterized by increased muscle bulk and strength. Affected individuals are exceptionally strong.
序列相似性	Belongs to the TGF-beta family.
细胞定位	Secreted.

图片



LANES:

1. MW marker – 14, 21, 31, 45, 66, 97 kDa
2. reduced and heated sample, 7µg/lane
3. non-reduced and non-heated sample, 7µg/lane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors