

Recombinant Human NELL1 protein ab116955

1 图像

描述

产品名称	重组人NELL1蛋白
表达系统	Wheat germ
Accession	<u>Q92832</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	

MPMDLILVVWFVCVCTARTVVGFGMDPDLQMDIVTELDLVNTT
LGVAQVSG
MHNASKAFLFQDIEREIHAAPHVSEKLIQLFRNKSEFTILAT
VQQKPSTS
GVILSIRELEHSYFELESSGLRDEIRYHYIHNGKPRTEALPY
RMADGQWH
KVALSVSASHLLLHVDCNRIYERVIDPPDTNLPPGINLWLGQ
RNQKHGLF
KGIIQDGKIIFMPNGYITQCPNLNHTCPTCSDFLSLVQGIMD
LQELLAKM
TAKLNYAETRLSQLENCHCEKTCQVSGLLYRDQDSWVDGDHC
RNCTCKSG
AVECRRMSCPPLNCSPDSL PVHIAGQCCKVCRPKCIYGGKVL
AEGQRILT
KSCRECRGGVLVKITEMCPPLNCSEKDHL PENQCCRVCRGH
NFCAEGPK
CGENSECKNWNTKATCECKSGYISVQGDSAYCEDIDECAAKM
HYCHANTV
CVNLPGLYRCDVPGYIRVDDFSCTEHDECGSGQHNCDENAI
CTNTVQGH
SCTCKPGYVGNGTICRAFCEEGRYGGTCVAPNKVCVPSGFT
GSHCEKDI
DECSEGIIECHNHSRCVNLPGWYHCECRSGFHDDGTYSLSGE
SCIDIDEC
ALRTHTCWNDSACINLAGGFDCLCPSGPSCSGDCPHEGGLKH
NGQVWTLK
EDRCSVCCKDGKIFCRRTACDCQNPSADLFCCPECDTRVTS
QCLDQNGH

KLYRSGDNWTHSCQQCRCLGEVDCWPLTCPNLSCEYTAILE
GECCPRCV
SDPCLADNITYDIRKTCLDSYGVSRLSGSVWTMAGSPCTTCK
CKNGRVCC SVDFECLQNN

预测分子量 115 kDa including tags
氨基酸 1 to 810

技术指标

Our **Abpromise guarantee** covers the use of **ab116955** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE
Western blot

形式 Liquid

制备和贮存

稳定性和存储 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.3% Glutathione, 0.79% Tris HCl


常规信息

功能 Involved in the control of cell growth and differentiation.

序列相似性 Contains 6 EGF-like domains.
Contains 1 TSP N-terminal (TSPN) domain.
Contains 5 VWFC domains.

细胞定位 Secreted.

图片



12.5% SDS-PAGE showing ab116955 at approximately 115.17kDa stained with Coomassie Blue.

SDS-PAGE - Recombinant Human NELL1 protein
(ab116955)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors