

Recombinant Human NHP2L1 protein ab101947

1 图像

描述

产品名称	重组人NHP2L1蛋白
纯度	> 95 % SDS-PAGE. ab101947 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>P55769</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMGSHMTEADVNPKAYPLADA HLTKKLLDLV QQSCNYKQLRKGANEATKTLNRGISEFIVMAADAEPLEIILH LPLLCEDK NVPYVFVRSKQALGRACGVS RPVIACSVTIKEGSQKQIQS IQSIERL LV
预测分子量	17 kDa including tags
氨基酸	1 to 128
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab101947** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

pH: 7.50


Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

相关性 Originally named because of its sequence similarity to the *Saccharomyces cerevisiae* NHP2 (non histone protein 2), this protein appears to be a highly conserved nuclear protein that is a component of the [U4/U6.U5] tri snRNP. It binds to the 5' stem loop of U4 snRNA. Two transcript variants encoding the same protein have been found for this gene.

细胞定位 Nucleus > nucleolus. Note: Concentrated in the dense fibrillar component of the nucleolus.

图片



15% SDS-PAGE analysis of 3µg ab101947.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors