

Recombinant Human NNT1 protein (denatured) ab111634

1 图像

描述	
产品名称	重组人NNT1蛋白(denatured)
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>Q9UBD9</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMLNRTGDPGPGPSIQKTYDLTR YLEHQLRS LAGTYLNYLGPPFNEPDFNPPRLGAETLPRATVDLEVWRSLN DKLRLTQN YEAYSHLLCYLRGLNRQAATAELRRSLAHFCTSLQGLLGSIA GVMAALGY PLPQPLPGTEPTWTPGPAHSDFLQKMDDFWLLKELQTLWRS AKDFNRLK KKMQPPAAAVTLHLGAHGF
预测分子量	25 kDa including tags
氨基酸	28 to 225
标签	His tag N-Terminus
描述	重组人NNT1蛋白

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab111634</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid

制备和贮存	
-------	--

## 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 3.50

Constituents: 0.59% Sodium citrate, 2.4% Urea, 10% Glycerol (glycerin, glycerine)

## 常规信息

### 功能

Cytokine with B-cell stimulating capability. Binds to and activates the ILST/gp130 receptor.

### 组织特异性

Expressed predominantly in lymph nodes, spleen, peripheral blood lymphocytes, bone marrow, and fetal liver.

### 疾病相关

Defects in CLCF1 are the cause of cold-induced sweating syndrome type 2 (CISS2) [MIM:610313]. Cold-induced sweating syndrome (CISS) is an autosomal recessive disorder characterized by profuse sweating induced by cool surroundings (temperatures of 7 to 18 degrees Celsius). Additional abnormalities include a high-arched palate, nasal voice, depressed nasal bridge, inability to fully extend the elbows and kyphoscoliosis.


### 序列相似性

Belongs to the IL-6 superfamily.

### 细胞定位

Secreted.

## 图片



15% SDS-PAGE showing ab111634 at approximately 24.6kDa (3μg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors