

Recombinant Human NUR77 protein ab114315

1 图像

描述

产品名称	重组人NUR77蛋白
表达系统	Wheat germ
Accession	<u>P22736</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<p>MPCIQAQYGTAPSPGPRDHLASDPLTPEFIKPTMDLASPEA APAAPTAL PSFSTFMDGYTGEFDTFLYQLPGTVQPCSSASSSASSTSSSS ATSPASAS FKFEDEFQVYGCYPGPLSGPVDEALSSSGSDYYGSPCSAPSPS TPSFQPPQ LSPWDGSFGHFSPSQTYEGLRAWTEQLPKASGPPQPPAFFSF SPPTGPSP SLAQSPKLKLFPSQATHQLGEGESYSMPATAFPGLAPTPHLEG SGILDTPV TSTKARSGAPGGSEGRCAVCGDNASCQHYGVRTCEGCKGFFK RTVQKNAK YICLANKDCPVDKRRRNRCQFCRFQKCLAVGMVKEVVRTDSL KGRRRLP SKPKQPPDASPANLLTSLVRAHLDSGPSTAKLDYSKFQELVL PHFGKEDA GDVQQFYDLLSGSLEVIRKWAEEKIPGFAELSPADQDLLLESA FLELFILR LAYRSKPGEGKLIFCSGLVLHRLQCARGFGDWIDSILAFSRS LHSLLVDV PAFACLSALVLTDRHGLQEPRRVEELQNRIASCLKEHVAAV AGEPQPAS CLSRLGKLPRLTCTQGLQRIFYLKLEDLVPPPPPIIDKIF MDTLPF</p>
预测分子量	92 kDa including tags
氨基酸	1 to 598

技术指标

Our **Abpromise guarantee** covers the use of **ab114315** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	SDS-PAGE
	Western blot
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

常规信息

功能	Orphan nuclear receptor. May act concomitantly with NURR1 in regulating the expression of delayed-early genes during liver regeneration. Binds the NGF1-B response element (NBRE) 5'-AAAAGGTCA-3' (By similarity). May inhibit NF-kappa-B transactivation of IL2.
组织特异性	Fetal muscle and adult liver, brain and thyroid.
序列相似性	Belongs to the nuclear hormone receptor family. NR4 subfamily. Contains 1 nuclear receptor DNA-binding domain.
细胞定位	Nucleus.

图片



SDS-PAGE analysis of ab114315 on a 12.5% gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors