

Recombinant Human NXT1 protein ab107147

1 图像

描述

| | |
|-----------|--|
| 产品名称 | 重组人NXT1蛋白 |
| 纯度 | > 85 % SDS-PAGE. Purified by using conventional chromatography |
| 表达系统 | Escherichia coli |
| Accession | <u>Q9UKK6</u> |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Human |
| 序列 | MGSSHHHHHSSGLVPRGSHMASVDFKTYVDQACRAAEFF VNVYYTTMDK RRRLLSRLYMGATLVWNGNAVSGQESLSEFFEMLPSEFQI SVVDCQPV HDEATPSQTTVLVVICGSVKFEGNKQRDFNQNFILTAQASPS NTVWKIAS DCFRFQDWAS |
| 预测分子量 | 18 kDa including tags |
| 氨基酸 | 1 to 160 |
| 标签 | His tag N-Terminus |

技术指标

Our **Abpromise guarantee** covers the use of **ab107147** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|-----|-------------------------------|
| 应用 | Mass Spectrometry SDS-PAGE |
| 质谱法 | MALDI-TOF |
| 形式 | Liquid |

制备和贮存

稳定性和存储 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

| | |
|--------------|---|
| 功能 | Stimulator of protein export for NES-containing proteins. Also plays a role in the nuclear export of U1 snRNA, tRNA, and mRNA. The NXF1-NXT1 heterodimer is involved in the export of HSP70 mRNA in conjunction with THOC4 and THOC5. |
| 序列相似性 | Contains 1 NTF2 domain. |
| 细胞定位 | Nucleus. Cytoplasm. Shuttles between the nucleus and the cytoplasm. |

图片



ab107147 analysed by SDS PAGE.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors