

Recombinant Human OTUB2 protein ab106896

1 图像

描述	
产品名称	重组人OTUB2蛋白
纯度	> 95 % SDS-PAGE. ab106896 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q96DC9</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMSETSFNLISEKCDILSILRDH PENRIYRR KIEELSKRFTAIRKTKGDGNCFYRALGYSYLESLLGKSREIF KFKERVLP TPNDLLAAGFEEHKFRNFFNAFYSVVELVEKDGSVSSLLKVF NDQSASDH IVQFLRLLTSAFIRNRADFFRHFIDEEMDIKDFCTHEVEPMA TECDHIQI TALSQALSIALQVEYVDEMDTALNHHVFPEAATPSVYLLYKT SHYNILYA ADKH
预测分子量	29 kDa including tags
氨基酸	1 to 234
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab106896 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

常规信息

功能

Hydrolase that can remove conjugated ubiquitin from proteins in vitro and may therefore play an important regulatory role at the level of protein turnover by preventing degradation. Mediates deubiquitination of both 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains, with a preference for 'Lys-63'-linked polyubiquitin chains.


组织特异性

Widely expressed. Expressed at higher level in brain.

序列相似性

Belongs to the peptidase C65 family.
Contains 1 OTU domain.

图片



3ug by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human OTUB2 protein
(ab106896)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors