

Recombinant Human Pallidin protein ab124579

1 图像

描述	
产品名称	重组人Pallidin蛋白
纯度	> 90 % SDS-PAGE. ab124579 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q9UL45</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MSVPGPSSPD GALTRPPYCL EAGEPTPGLS DTSPDEGLIE DLTIEDKAVE QLAEGLLSHY LPDLQRSKQA LQELTQNQVV LLDLTLEQEIS KFKECHSMLD INALFAEAKH YHAKLVNIRK EMLMLHEKTS KLKKRALKLQ QKRQKEELER EQQREKEFER EKQLTARPAK RM
预测分子量	22 kDa including tags
氨基酸	1 to 172
标签	His tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab124579 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

功能

Involved in the development of lysosome-related organelles, such as melanosomes and platelet-dense granules. May play a role in intracellular vesicle trafficking, particularly in the vesicle-docking and fusion process.

组织特异性

Widely expressed.


翻译后修饰

Phosphorylated.

细胞定位

Cytoplasm. Endomembrane system. It can exist as a soluble protein as well as a peripheral membrane protein.

图片



15% SDS-PAGE showing ab124579 (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors