

Product datasheet

Recombinant Human PAR1/Thrombin Receptor protein ab114346

1 图像

描述

产品名称	重组人PAR1/Thrombin Receptor蛋白
表达系统	Wheat germ
Accession	P25116
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	SFLLRNPNDKYEPFWEDEEKNESGLTEYRLVSIKSSPLQKQ LPAFISED ASGYLTSSWLT
预测分子量	32 kDa including tags
氨基酸	42 to 102

技术指标

Our **Abpromise guarantee** covers the use of **ab114346** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Western blot ELISA
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

常规信息

功能	High affinity receptor for activated thrombin coupled to G proteins that stimulate phosphoinositide hydrolysis. May play a role in platelets activation and in vascular development.
组织特异性	Platelets and vascular endothelial cells.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
翻译后修饰	A proteolytic cleavage generates a new N-terminus that functions as a tethered ligand. Phosphorylated; probably mediating desensitization prior to the uncoupling and internalization of the receptor.
细胞定位	Cell membrane.

图片



ab114346 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors