

Product datasheet

Recombinant Human PBLD protein ab101217

1 图像

描述

产品名称	重组人PBLD蛋白
纯度	> 95 % SDS-PAGE. ab101217 was purified using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
表达系统	Escherichia coli
Accession	<u>P30039</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<p>MGSSHHHHHSSGLVPRGSHMKLPIFIADAFTARAFRGNP AAVCLLENEL DEDMHQKIAREMNLSETAFIRKLHPTDNFAQSSCFGLRWFTP ASEVPLCG HATLASAAVL FHKIKNMNSTLTFVTLSGELRARRAEDGIVLD LPLYPAHP QDFHEVEDLIKTAIGNTLVQDICYSPDTQKLLVRLSDVYNRS FLENLKVN TENLLQVENTGKVKGLILTLKGEPGGQTQAFDFYSRYFAPWV GVAEDPVT GSAHAVLSSYSQHLGKKEMHAFQCShRGGELGISLRPDGRV DIRGGAAV VLEGLTLTA</p>
预测分子量	34 kDa
氨基酸	1 to 288
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab101217** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	Mass Spectrometry

质谱法 MALDI-TOF
形式 Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

序列相似性

Belongs to the phzF family.

图片



ab101217 at 3 µg analysed by 15% SDS PAGE.

SDS-PAGE - Recombinant Human PBLD protein
(ab101217)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors