

Recombinant human PDE4B protein ab125582

5 图像

描述

产品名称	重组人PDE4B蛋白
生物活性	The specific activity of ab125582 was determined to be 538 nmol/min/mg.
纯度	>= 70 % Densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<u>Q07343</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	78 kDa including tags
氨基酸	305 to 736
标签	GST tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab125582** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.88% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---


常规信息

功能	Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes. May be involved in mediating central nervous system effects of therapeutic agents ranging from antidepressants to antiasthmatic and anti-inflammatory agents.
组织特异性	Expressed in brain, heart, lung and skeletal muscle.
通路	Purine metabolism; 3',5'-cyclic AMP degradation; AMP from 3',5'-cyclic AMP: step 1/1.
序列相似性	Belongs to the cyclic nucleotide phosphodiesterase family. PDE4 subfamily.


图片




The specific activity of PDE4B (ab125582) was determined to be 470 nmol/min/mg as per activity assay protocol




SDS PAGE analysis of ab125582



SDS PAGE analysis of ab125582



SDS Page analysis of ab125582



The specific activity of ab125582 was determined to be 538 nmol/min/mg.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors