

Recombinant Human PFD6 protein ab104644

1 图像

描述	
产品名称	重组人PFD6蛋白
纯度	> 95 % SDS-PAGE. ab104644 is purified the proteins by using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
表达系统	Escherichia coli
Accession	<u><b>O15212</b></u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<b>MGSSHHHHHHSSGLVPRGSHMAELIQKKLQGEVEKYQQLQ KDLSKSMGR QKLEAQLTENNIVKEELALLDGSNVVFKLLGPVLVKQELGEA RATVGKRL DYITAEIKRYESQLRDLERQSEQQRETLAQLQQEFQRAQAAK AGAPGKA</b>
预测分子量	17 kDa including tags
氨基酸	1 to 129
标签	His tag N-Terminus
技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab104644</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	

## 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 常规信息


### 功能

Binds specifically to cytosolic chaperonin (c-CPN) and transfers target proteins to it. Binds to nascent polypeptide chain and promotes folding in an environment in which there are many competing pathways for nonnative proteins.

### 序列相似性

Belongs to the prefoldin subunit beta family.

## 图片



15% SDS-PAGE showing ab104644 at approximately 16.7kDa (3μg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors