

# Recombinant human PI3K (p110 delta/p85 alpha) protein (Active) ab125633

8 References 2 图像

描述	
产品名称	重组人PI3K (p110 delta/p85 alpha)蛋白(Active)
生物活性	The specific activity of ab125633 was equivalent to 2,514 nmol/min/mg as per activity assay protocol.
纯度	> 75 % Densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<u><b>O00329</b></u> <u><b>P27986</b></u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
氨基酸序列 1	
种属	Human
预测分子量	111 kDa
氨基酸	1 to 1044
标签	His tag N-Terminus
额外的序列信息	NM_005026.
氨基酸序列 2	
种属	Human
预测分子量	86 kDa
氨基酸	2 to 724
标签	His tag N-Terminus
额外的序列信息	NM_181253.
技术指标	

Our **Abpromise guarantee** covers the use of **ab125633** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot Functional Studies SDS-PAGE
形式	Liquid
补充说明	Kinase assay.


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.002% PMSF, 0.82% Sodium phosphate, 0.004% DTT, 25% Glycerol, 1.74% Sodium chloride  This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

图片



The specific activity of ab125633 was equivalent to 2,514 nmol/min/mg as per activity assay protocol.



SDS-PAGE - Recombinant human PI3K (p110 delta/p85 alpha) protein (Active) (ab125633)

Purity of ab125633 was determined to be >75% by densitometry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors