

Recombinant Human pNO40 protein ab104634

1 图像

描述	
产品名称	重组人pNO40蛋白
纯度	> 90 % SDS-PAGE. ab104634 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q9NP64</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMGSMNSGRPETMENLPALYT IFQGEVAMVT DYGAFIKIPGCRKQGLVHRTHMSSCRVDKPSEIVDVGDKVWV KLIGREMK NDRIKVSLSMKVVNQGTGKDLDPNNVIEQEERRRRSFQDYT GQKITLEA VLNTTCKKCGCKGHFAKDCFMQPGGTKYSLIPDEEEEKEEAK SAEFKPD PTRNPSRKRKKEKKKKKHRDRKSSSDSDSSDESDTGKRARHT SKDSKAAK KKKKKKKHKKKHKE
预测分子量	30 kDa including tags
氨基酸	1 to 241
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab104634 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid

补充说明

Previously labelled as ZCCHC17.

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00174% PMSF, 0.0308% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

序列相似性


Contains 1 CCHC-type zinc finger.

Contains 1 S1 motif domain.

细胞定位

Nucleus, nucleolus.

图片



15% SDS-PAGE showing ab104634 (3 µg).

Note: Molecular weight on SDS-PAGE will appear higher than the predicted molecular weight of approximately 30 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors