

Recombinant human PPIH protein ab78874

1 图像

描述

产品名称	重组人PPIH蛋白
生物活性	Specific activity is > 220nmol/min/mg, and is defined as the amount of enzyme that cleaves 1umole of suc-AAPF-pNA per minute at 37°C in Tris-Hcl pH8.0 using chymotrypsin.
纯度	> 95 % SDS-PAGE. ab78874 is purified using conventional chromatography techniques.
内毒素水平	< 1.000 Eu/μg
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAVANSSPVNPVVFFDVSIGGQEVGRMKIELFADVVPKTAEN FRQFCTGE FRKDGVPIGYKGSTFHRVIKDFMIQGGDFVNGDGTGVASIYR GPFADENF KLRHSAPGLLSMANS GPSTNGCQFFITCSKCDWLDGKHVVF G KIIDGLLV MRKIENVPTGPNNKPKLPVVISQCGEM

技术指标

Our **Abpromise guarantee** covers the use of **ab78874** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
---------------	--

pH: 7.40


Constituents: PBS, 10% Glycerol (glycerin, glycerine)

This product is an active protein and may elicit a biological response in vivo, handle with caution.

常规信息

功能	PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. Participates in pre-mRNA splicing. May play a role in the assembly of the U4/U5/U6 tri-snRNP complex. May act as a chaperone.
序列相似性	Belongs to the cyclophilin-type PPlase family. PPlase H subfamily. Contains 1 PPlase cyclophilin-type domain.
细胞定位	Nucleus speckle. Cytoplasm. Colocalizes with spliceosomal snRNPs. A small proportion may also be cytoplasmic.

图片



15% SDS-PAGE showing ab78874 at approximately 19kDa (3 μ g).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors