

Recombinant Human Psoriasin protein ab105577

1 References 2 图像

描述	
产品名称	重组人Psoriasin蛋白
纯度	> 90 % SDS-PAGE. ab105577 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	P31151
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMSNTQAERSIIGMIDMFHKY TRRDDKIDKP SLLTMMKENFPNFLSACDKKGTNYLADVFEKKDKNEDKKIDF SEFLSLLG DIATDYHKQSHGAAPCSGGSQ
预测分子量	14 kDa including tags
氨基酸	1 to 101
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab105577 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

组织特异性

Fetal ear, skin, and tongue and human cell lines. Highly up-regulated in psoriatic epidermis. Also highly expressed in the urine of bladder squamous cell carcinoma (SCC) bearing patients.


序列相似性

Belongs to the S-101 family.
Contains 2 EF-hand domains.

细胞定位


Cytoplasm. Secreted. Secreted by a non-classical secretory pathway.

图片



SDS-PAGE - Recombinant Human Psoriasin protein (ab105577)

15% SDS-PAGE analysis of ab105577 (3ug)



Western blot - Recombinant Human Psoriasin protein (ab105577)

Anti-Psoriasin antibody ([ab83534](#)) at 1 µg/ml + Recombinant Human Psoriasin protein (ab105577) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ([ab97080](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 1 minute

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors