

Recombinant Human Rab24 protein ab113161

1 图像

描述	
产品名称	重组人Rab24蛋白
纯度	> 90 % SDS-PAGE. ab113161 was purified by using conventional chromatography techniques
表达系统	Escherichia coli
Accession	<u>Q969Q5</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMSGQRVDVKVVMLGKE YVGKTSLVER YVHDRFLVGPYQNTIGAAFVAKVMSVGDRVTTLGIWDTAGSE RYEAMSRI YYRGAKAAIVCYDLTDSSSFERAKFWVKELRSLEEGCQIYLC GTKSDLLE EDRRRRRVDFHDVQDYADNIKAQLFETSSKTGQSVDEL FQKV AEDYVSVA AFQVMTEDKGVDLGQKPNPYFYSCCHH
预测分子量	26 kDa including tags
氨基酸	1 to 227
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab113161 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	

稳定性和存储	<p>Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.</p> <p>pH: 8.00</p> <p>Constituents: 0.03% DTT, 0.32% Tris HCl, 0.03% EDTA, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride</p>
---------------	---

常规信息

功能	May be involved in autophagy-related processes.
序列相似性	Belongs to the small GTPase superfamily. Rab family.
翻译后修饰	Isoprenylation is inefficient compared to other Rab family members.
细胞定位	Cytoplasm > cytosol. Membrane. Only about 20-25% is recovered in the particulate fraction.

图片



ab113161 at 3 µg analysed by 15% SDS PAGE.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors