

### Recombinant human Rab9 protein ab90106

#### 描述

---

产品名称	重组人Rab9蛋白
生物活性	Activity: 100 pmol of protein can bind > 80 pmol of GDP.
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human

#### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab90106** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
形式	Liquid

#### 制备和贮存

---

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.00044% GDP, 0.095% Magnesium chloride, 0.0156% Beta mercaptoethanol, 0.79% Tris HCl, 0.58% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---

#### 常规信息

---

功能	Involved in the transport of proteins between the endosomes and the trans Golgi network.
序列相似性	Belongs to the small GTPase superfamily. Rab family.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors