

Recombinant Human RCAS1 protein ab79243

1 图像

描述

产品名称	重组人RCAS1蛋白
纯度	> 90 % SDS-PAGE. ab79243 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> MGSSHHHHHH SSGLVPRGSH MRSGRGRKLS GDQITLPTTV DYSSVPKQTD VEEWTSWDED APTSVKIEGG NGNVATQQNS LEQLEPDYFK DMTPTIRKTQ KIVIKKREPL NFGIPDGSTG FSSRLAATQD LPFIHQSSSEL GDLDTWQENT NAWEEEEEDAA WQAEVLRQQ KLADREKRAA EQQRKKMEKE AQRLMKKEQN KIGVKLS </pre>

技术指标

Our **Abpromise guarantee** covers the use of **ab79243** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	<p>Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.</p> <p>pH: 8.00</p> <p>Constituents: 0.0154% DTT, 0.242% Tris, 0.0584% EDTA, 10% Glycerol (glycerin, glycerine)</p>
--------	--

常规信息

功能	May participate in suppression of cell proliferation and induces apoptotic cell death through activation of interleukin-1-beta converting enzyme (ICE)-like proteases.
组织特异性	Widely expressed. Expressed in ovary, testis, prostate, thymus, muscle and heart, but not in small intestine, colon, lymph nodes, or peripheral blood lymphocytes. The protein is not detected in any of the above organs.
结构域	The coiled coil domain is necessary for the homodimerization.
细胞定位	Golgi apparatus membrane. According to PubMed:10426319, it also exists as a soluble form which has the same biological activities. The existence of such soluble form is however uncertain.

图片



15% SDS-PAGE showing ab79243 (3 μ g).

Note: Real molecular weight on SDS-PAGE is shifted up compared to predicted molecular weight.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors