

Recombinant Human RTN4IP1 protein ab116191

1 图像

描述	
产品名称	重组人RTN4IP1蛋白
纯度	> 95 % SDS-PAGE. ab116191 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	<u>Q8WWV3</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHH SSGLVP RGSHMGS TVMPAWVIDKYGKNEVLRFTQNNMMPIIHYPNEVIVKV HAASVNPIDVNM RSGYGATALNMKRDPLHVKIKGEEFPLTLGRDVSGVVM ECGLDVKYFKPG DEVWAAVPPWKQGTLSEFVVVSGNEVSHKPKSLTHTQA ASLPYVALTAWS AINKVGGLNDKNCTGKRVLILGASGGVGTFAIQVMKAW DAHVTAVCSQDA SELVRKLGADDVIDYKSGSVEEQLKSLKPFDFILDNVG GSTETWAPDFLK KWSGATYVTLVTPFLLNMDRLGIADGMLQTGVTVGSKA LKHFWKGVHYRW AFFMASGPCLDDIAELVDAGKIRPVIEQTFPFSKVPEA FLKVERGHARGK TVINVV
预测分子量	41 kDa including tags
氨基酸	41 to 396
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab116191** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry
	SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine)
--------	---

常规信息

功能	Appears to be a potent inhibitor of regeneration following spinal cord injury.
组织特异性	Widely expressed in mitochondria-enriched tissues. Found in heart, muscle, kidney, liver, brain and placenta.
序列相似性	Belongs to the zinc-containing alcohol dehydrogenase family. Quinone oxidoreductase subfamily.
细胞定位	Mitochondrion.

图片



15% SDS-PAGE analysis of 3 µg ab116191.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors