

# Recombinant Human RUNX3 protein ab95893

## 1 图像

### 描述

产品名称	重组人RUNX3蛋白
纯度	> 95 % SDS-PAGE. ab95893 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u><a href="#">Q13761</a></u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MRSMVDVLAD HAGELVRTDS PNFLCSVLPV HWRCNKTLPV AFKVVALGDV PDGTVVTVMA GNDENYSAEL RNASAVMKNQ VARFNDLRFV GRSGRGKSFT LTITVFTNPT QVATYHRAIK VTVDGPREPR RHRQK
预测分子量	17 kDa including tags
氨基酸	53 to 186
标签	His tag N-Terminus

### 技术指标

Our **Abpromise guarantee** covers the use of **ab95893** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF-TOF
形式	Liquid

### 制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
--------	--


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 常规信息

<b>功能</b>	CBF binds to the core site, 5'-PYGPYGGT-3', of a number of enhancers and promoters, including murine leukemia virus, polyomavirus enhancer, T-cell receptor enhancers, Ick, IL-3 and GM-CSF promoters.
<b>序列相似性</b>	Contains 1 Runt domain.
<b>结构域</b>	A proline/serine/threonine rich region at the C-terminus is necessary for transcriptional activation of target genes.
<b>翻译后修饰</b>	Phosphorylated on tyrosine residues by SRC. Phosphorylated by LCK and FYN.
<b>细胞定位</b>	Nucleus. Cytoplasm. The tyrosine phosphorylated form localizes to the cytoplasm.

## 图片



15% SDS-PAGE showing ab95893 at approximately 17.1kDa (3 $\mu$ g).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors