

Recombinant Human RXRG protein ab114438

1 图像

描述		
产品名称	重组人RXRG蛋白	
表达系统	Wheat germ	
Accession	<u>P48443</u>	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MYGNYSHFMKFPAGYGGSPGHTGSTSMSPSAALSTGKPMDSH PSYTDTPV SAPRTL SAVGTPLNALGSPYRVITSAMGPPSGALAAPPGINL VAPPSSQL NVVNSVSSSEDIKPLPGLPGIGNMYPSTSPGSLVKHICAIC GDRSSGKH YGVYSCEGCKGFFKRTIRKDLIYTCDNKDCLIDKRQRNRCQ YCRYQKCL VMGMKREAVQEERQRSRERAESAEACATSGHEDMPVERILEA ELAVEPKT ESYGDMNMENSTNDPVTNICHAAADKQLFTLVWAKRIPHFSD LTLEDQVI LLRAGWNELLIASFSHRSVSVQDGILLATGLHVHRSSAHSAG VGSIFDRV LTELVSKMKDMQMDKSELGCLRAIVLFNPDAKGLSNPSEVET LREKVYAT LEAYTKQKYPEQGRFAKLLLRPALRSIGLKCLEHLFFFKL IGDTPIDT FLMEMLETPLQIT	
预测分子量	77 kDa including tags	
氨基酸	1 to 463	

技术指标


Our **Abpromise guarantee** covers the use of **ab114438** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
----	-------

	Western blot
	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as Retinoid X Receptor gamma.
制备和贮存	
稳定性和存储	<p>Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.</p> <p>pH: 8.00</p> <p>Constituents: 0.3% Glutathione, 0.79% Tris HCl</p>
常规信息	
功能	Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes. The RAR/RXR heterodimers bind to the retinoic acid response elements (RARE) composed of tandem 5'-AGGTCA-3' sites known as DR1-DR5. The high affinity ligand for RXRs is 9-cis retinoic acid.
序列相似性	<p>Belongs to the nuclear hormone receptor family. NR2 subfamily.</p> <p>Contains 1 nuclear receptor DNA-binding domain.</p>
结构域	Composed of three domains: a modulating N-terminal domain, a DNA-binding domain and a C-terminal ligand-binding domain.
细胞定位	Nucleus.

图片



ab114438 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors