

Recombinant Human S100A11 protein ab55569

1 图像

描述

产品名称	重组人S100A11蛋白
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human

技术指标

Our **Abpromise guarantee** covers the use of **ab55569** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. pH: 7.40 Constituents: PBS, 50% Glycerol
--------	---

常规信息

功能	Facilitates the differentiation and the cornification of keratinocytes.
序列相似性	Belongs to the S-100 family. Contains 2 EF-hand domains.
翻译后修饰	Phosphorylation at Thr-10 by PRKCA significantly suppresses homodimerization and promotes association with NCL/nucleolin which induces nuclear translocation.
细胞定位	Cytoplasm. Nucleus.

图片



ab55569 on SDS-PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors