

Recombinant Human S100P protein ab63309

1 图像

描述	
产品名称	重组人S100P蛋白
纯度	> 95 % SDS-PAGE. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.
表达系统	Escherichia coli
Accession	<u>P25815</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MTELEAAMGM IIDVFSRYSG SEGSTQTLTK GELKVLMEKE LPGFLQSGKD KDAVDKLLKD LDANGDAQVD FSEFIVFVAA ITSACHKYFE KAGLK
预测分子量	12 kDa including tags
氨基酸	1 to 95
标签	His tag N-Terminus

技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab63309 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot ELISA SDS-PAGE
形式	Lyophilized
补充说明	ab63309 has an N-terminal His-tag (10 addiitonal amino acids).

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.

Constituent: 0.134% PBS

复溶

Add deionized water to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely.

常规信息

组织特异性

Up-regulated in various pancreatic ductal adenocarcinomas and pancreatic intraepithelial neoplasias.


序列相似性

Belongs to the S-100 family.
Contains 2 EF-hand domains.

细胞定位

Nucleus. Cytoplasm. Colocalizes with S100PBP in the nucleus.

图片



1. MW marker – 14, 21, 31, 45, 66, 97 kDa
2. reduced and heated sample, 5µg/lane
3. non-reduced and non-heated sample, 5µg/lane

SDS-PAGE - Recombinant Human S100P protein
(ab63309)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors