

Recombinant Human SAE1 protein ab96772

1 References 1 图像

描述	
产品名称	重组人SAE1蛋白
纯度	> 90 % SDS-PAGE. ab96772 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q9UBE0</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MHHHHHHMAS MTGGQQMGRD LYDDDDKDRW GSMVEKEEAG GGISEEEAAQ YDRQIRLWGL EAQKRLRASR VLLVGLKGLG AEIAKNLILA GVKGLTMLDH EQVTPEDPGA QFLIRTGSVG RNRAEASLER AQNLNPMVDV KVDTEIEKK PESFFTQFDA VCLTCCSRDV IVKVDQICHK NSIKFFTGDV FGYHGYTFAN LGEHEFVEEK TKVAKVSQGV EDGPDTKRAK LDSSETTMVK KKVVFPCVKE ALEVDWSSEK AKAALKRTTS DYFLLQVLLK FRTDKGRDPS SDTYEEDSEL LLQIRNDVLD SLGISPDLLP EDFVRYCFSE MAPVCAVVG ILAQEIVKAL SQRDPPHNNF FFFDGMKGNG IVECLGPK
预测分子量	42 kDa including tags
氨基酸	1 to 346
标签	His tag N-Terminus , T7 tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab96772</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE

质谱法	MALDI-TOF-TOF
形式	Liquid
<b>制备和贮存</b>	
稳定性和存储	<p>Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.</p> <p>pH: 8.00</p> <p>Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)</p>
<b>常规信息</b>	
功能	The heterodimer acts as a E1 ligase for SUMO1, SUMO2, SUMO3, and probably SUMO4. It mediates ATP-dependent activation of SUMO proteins followed by formation of a thioester bond between a SUMO protein and a conserved active site cysteine residue on UBA2/SAE2.
组织特异性	Expression level increases during S phase and drops in G2 phase (at protein level).
通路	Protein modification; protein sumoylation.
序列相似性	Belongs to the ubiquitin-activating E1 family.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Nucleus.

## 图片



15% SDS-PAGE showing ab96772 at approximately 42.2kDa (3μg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors