

Recombinant Human Serpin A5 protein ab152068

描述	
产品名称	重组人Serpin A5蛋白
纯度	> 95 % SDS-PAGE. The purity of ab152068 is greater than 95%, as determined by SEC-HPLC and reducing SDS-PAGE. It is lyophilized from a 0.2 μM filtered solution.
内毒素水平	< 1.000 Eu/μg
表达系统	HEK 293 cells
Accession	<u>P05154</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	HRHHPREMKKRVEDLHVGATVAPSSRRDFTFDLYRALASAAP SQNIFFSP VSISMSLAMLSLGAGSSTKMQILEGLGLNLQKSSEKELHRGF QQLLQELN QPRDGFQLSLGNALFTDLVVDLQDTFVSAMKTLYLADTFPTN FRDSAGAM KQINDYVAKQTKGKIVDLLKNLDSNAVVMVNYIFFKAKWET SFNHKGTQ EQDFYVTSETVVRVPMMSREDQYHYLLDRNLSCRVVGVPYQG NATALFIL PSEGKMQQVENGLSEKTLRKWLKMFKKRQLELYLPKFSIEGS YQLEKVLP SLGISNVFTSHADLSGISNHSNIQVSEMVHKAVVEVDESGTR AAAATGTI FTFRSARLNSQRLVFNRPFILMFIVDNNILFLGKVNRPVDHHH HHH
预测分子量	45 kDa including tags
氨基酸	20 to 406
标签	His tag C-Terminus
技术指标	

Our **Abpromise guarantee** covers the use of **ab152068** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE HPLC
形式	Lyophilized
补充说明	Previously labelled as Protein C inhibitor.

制备和贮存

稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term. pH: 5.50 Constituents: 0.39% MES, 0.88% Sodium chloride
复溶	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in 1X PBS. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

常规信息

功能	Inhibits activated protein C as well as plasminogen activators.
组织特异性	Expressed by the liver and secreted in plasma.
序列相似性	Belongs to the serpin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors